

For Research Use Only

# INHBC Polyclonal antibody, PBS Only

Catalog Number: 24196-1-PBS



## Basic Information

Catalog Number:

24196-1-PBS

Size:

100ug, Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21358

GenBank Accession Number:

BC130326

GeneID (NCBI):

3626

UNIPROT ID:

P55103

Full Name:

inhibin, beta C

Calculated MW:

352 aa, 38 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data