

For Research Use Only

IL-28B/ IFNL3 Polyclonal antibody, PBS Only

Catalog Number: 24199-1-PBS



Basic Information

Catalog Number:

24199-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21491

GenBank Accession Number:

BC130316

GeneID (NCBI):

282617

UNIPROT ID:

Q8IZI9

Full Name:

interleukin 28B (interferon, lambda 3)

Calculated MW:

196 aa, 22 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using 24199-1-AP (IL28B/ IFNL3 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 24199-1-PBS in a different storage buffer formulation.