

For Research Use Only

# FOXB1 Polyclonal antibody

Catalog Number: 24285-1-AP



## Basic Information

Catalog Number:

24285-1-AP

Size:

150ul , Concentration: 900 µg/ml by Nanodrop and 513 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18079

GenBank Accession Number:

BC113710

GeneID (NCBI):

27023

UNIPROT ID:

Q99853

Full Name:

forkhead box B1

Calculated MW:

325 aa, 35 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HepG2 cells, mouse liver tissue

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 24285-1-AP (FOX B1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



mouse liver tissue were subjected to SDS PAGE followed by western blot with 24285-1-AP (FOX B1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.