

For Research Use Only

FOXB1 Polyclonal antibody

Catalog Number: 24285-1-AP



Basic Information

Catalog Number: 24285-1-AP	GenBank Accession Number: BC113710	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 513 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 27023	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: forkhead box B1	
Isotype: IgG	Calculated MW: 325 aa, 35 kDa	
Immunogen Catalog Number: AG18079	Observed MW: 35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells, mouse liver tissue
Species Specificity: human, mouse	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 24285-1-AP (FOX B1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



mouse liver tissue were subjected to SDS PAGE followed by western blot with 24285-1-AP (FOX B1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.