

For Research Use Only

LARGE Polyclonal antibody

Catalog Number: 24307-1-AP



Basic Information

Catalog Number:

24307-1-AP

Size:

150ul , Concentration: 300 ug/ml by Nanodrop and 233 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19458

GenBank Accession Number:

BC126404

GeneID (NCBI):

9215

UNIPROT ID:

O95461

Full Name:

like-glycosyltransferase

Calculated MW:

756 aa, 88 kDa

Observed MW:

65-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse kidney tissue, mouse brain tissue, mouse skeletal muscle tissue, NIH/3T3 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 24307-1-AP (LARGE Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.