

For Research Use Only

LARGE Polyclonal antibody

Catalog Number: 24307-1-AP



Basic Information

Catalog Number: 24307-1-AP	GenBank Accession Number: BC126404	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop and 233 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 9215	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: Like-glycosyltransferase	
Isotype: IgG	Calculated MW: 756 aa, 88 kDa	
Immunogen Catalog Number: AG19458	Observed MW: 65-70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: human, mouse	WB : mouse kidney tissue, mouse brain tissue, mouse skeletal muscle tissue, NIH/3T3 cells

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 24307-1-AP (LARGE Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.