

For Research Use Only

UGT2A1 Polyclonal antibody

Catalog Number: 24343-1-AP



Basic Information

Catalog Number: 24343-1-AP	GenBank Accession Number: BC157012	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop;	GeneID (NCBI): 10941	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: UDP glucuronosyltransferase 2 family, polypeptide A1	
Isotype: IgG	Calculated MW: 527 aa, 60 kDa	
Immunogen Catalog Number: AG19421	Observed MW: 50-53 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : SH-SY5Y cells, C6 cells
Species Specificity: human, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24343-1-AP (UGT2A1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.