

For Research Use Only

UGT2A1 Polyclonal antibody

Catalog Number: 24343-1-AP



Basic Information

Catalog Number:

24343-1-AP

Size:

150ul , Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19421

GenBank Accession Number:

BC157012

GeneID (NCBI):

10941

UNIPROT ID:

PODTE4

Full Name:

UDP glucuronosyltransferase 2 family, polypeptide A1

Calculated MW:

527 aa, 60 kDa

Observed MW:

50-53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : SH-SY5Y cells, C6 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24343-1-AP (UGT2A1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.