

For Research Use Only

DMRTA1 Polyclonal antibody

Catalog Number: 24356-1-AP



Basic Information

Catalog Number:

24356-1-AP

Size:

150ul , Concentration: 650 ug/ml by Nanodrop and 653 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18796

GenBank Accession Number:

BC130437

GeneID (NCBI):

63951

UNIPROT ID:

Q5VZB9

Full Name:

DMRT-like family A1

Calculated MW:

504 aa, 53 kDa

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

WB, IP, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse liver tissue,

IP : mouse liver tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 24356-1-AP (DMRTA1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.

IP result of anti-DMRTA1 (IP:24356-1-AP, 3ug; Detection:24356-1-AP 1:500) with mouse liver tissue lysate 3000ug.