

For Research Use Only

DMRTA1 Polyclonal antibody, PBS Only

Catalog Number: 24356-1-PBS



Basic Information

Catalog Number:

24356-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18796

GenBank Accession Number:

BC130437

GeneID (NCBI):

63951

UNIPROT ID:

Q5VZB9

Full Name:

DMRT-like family A1

Calculated MW:

504 aa, 53 kDa

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 24356-1-AP (DMRTA1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24356-1-PBS in a different storage buffer formulation.



IP result of anti-DMRTA1 (IP:24356-1-AP, 3ug; Detection:24356-1-AP 1:500) with mouse liver tissue lysate 3000ug. This data was developed using the same antibody clone with 24356-1-PBS in a different storage buffer formulation.