

For Research Use Only

ACSM1 Polyclonal antibody

Catalog Number: 24407-1-AP



Basic Information

Catalog Number:

24407-1-AP

Size:

150ul , Concentration: 310 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18475

GenBank Accession Number:

BC125178

GeneID (NCBI):

116285

UNIPROT ID:

Q08AH1

Full Name:

acyl-CoA synthetase medium-chain family member 1

Calculated MW:

577 aa, 65 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : K-562 cells, LNCaP cells, PC-3 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24407-1-AP (ACSM1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.