

For Research Use Only

Collagen Type IV Polyclonal antibody



Catalog Number: 24460-1-AP

1 Publications

Basic Information

Catalog Number:

24460-1-AP

Size:

150ul, Concentration: 360 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21443

GenBank Accession Number:

BC151220

GeneID (NCBI):

1282

UNIPROT ID:

P02462

Full Name:

collagen, type IV, alpha 1

Calculated MW:

1669 aa, 161 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IHC

Species Specificity:

human, rat, mouse

Cited Species:

rat

Notable Publications

Author	Pubmed ID	Journal	Application
Ping Yang	31572528	Exp Ther Med	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

