

For Research Use Only

Collagen Type IV Polyclonal antibody



Catalog Number: 24460-1-AP

1 Publications

Basic Information

Catalog Number: 24460-1-AP	GenBank Accession Number: BC151220	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 360 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 1282	
Source: Rabbit	Full Name: collagen, type IV, alpha 1	
Isotype: IgG	Calculated MW: 1669 aa, 161 kDa	
Immunogen Catalog Number: AG21443		

Applications

Tested Applications:
ELISA

Cited Applications:
IHC

Species Specificity:
human, rat, mouse

Cited Species:
rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ping Yang	31572528	Exp Ther Med	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data