

For Research Use Only

FBXL17 Polyclonal antibody, PBS Only

Catalog Number: 24476-1-PBS



Basic Information

Catalog Number:

24476-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20436

GenBank Accession Number:

BC126144

GeneID (NCBI):

64839

UNIPROT ID:

Q9UF56

Full Name:

F-box and leucine-rich repeat protein 17

Calculated MW:

701 aa, 76 kDa

Observed MW:

76 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

FBXL17 (F-box and leucine-rich repeat protein 17) is a substrate-recognition component of the SCF-type E3 ubiquitin ligase complex (SKP1-CUL1-F-box).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24476-1-AP (FBXL17 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24476-1-PBS in a different storage buffer formulation.