

For Research Use Only

LCK Polyclonal antibody, PBS Only

Catalog Number: 24497-1-PBS



Basic Information

Catalog Number:

24497-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21475

GenBank Accession Number:

BC013200

GeneID (NCBI):

3932

UNIPROT ID:

P06239

Full Name:

Lymphocyte-specific protein tyrosine kinase

Calculated MW:

539 aa, 56 kDa

Observed MW:

50-60 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 24497-1-AP (LCK Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24497-1-PBS in a different storage buffer formulation.