

For Research Use Only

# PPP2R3B Polyclonal antibody

Catalog Number: 24502-1-AP



## Basic Information

|   |  |  |
|---|--|--|
| <b>Catalog Number:</b><br>24502-1-AP  | <b>GenBank Accession Number:</b><br>BC009032   | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul , Concentration: 1000 µg/ml by Nanodrop and 547 µg/ml by Bradford method using BSA as the standard; | <b>GeneID (NCBI):</b><br>28227   | <b>Recommended Dilutions:</b><br>WB 1:500-1:2000             |
| <b>Source:</b><br>Rabbit  | <b>Full Name:</b><br>protein phosphatase 2 (formerly 2A), regulatory subunit B <sup>β</sup> , beta |  |
| <b>Isotype:</b><br>IgG  | <b>Calculated MW:</b><br>65 kDa  |  |
| <b>Immunogen Catalog Number:</b><br>AG21366   | <b>Observed MW:</b><br>65 kDa  |  |

## Applications

|  |   |
|--|---|
| <b>Tested Applications:</b><br>WB, ELISA | <b>Positive Controls:</b><br>WB : HL-60 cells, HeLa cells |
| <b>Species Specificity:</b><br>human     |   |

## Background Information

### Storage


**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HL-60 cells were subjected to SDS PAGE followed by western blot with 24502-1-AP (PPP2R3B Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.