

For Research Use Only

NEU2 Polyclonal antibody, PBS Only

Catalog Number: 24523-1-PBS



Basic Information

Catalog Number:

24523-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17386

GenBank Accession Number:

BC107053

GeneID (NCBI):

4759

UNIPROT ID:

Q9Y3R4

Full Name:

sialidase 2 (cytosolic sialidase)

Calculated MW:

380 aa, 42 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-12 cells were subjected to SDS PAGE followed by western blot with 24523-1-AP (NEU2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24523-1-PBS in a different storage buffer formulation.