

For Research Use Only

IL-25 Polyclonal antibody

Catalog Number: 24538-1-AP

1 Publications



Basic Information

Catalog Number:

24538-1-AP

Size:

150ul, Concentration: 187 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21441

GenBank Accession Number:

BC069565

GeneID (NCBI):

64806

UNIPROT ID:

Q9H293

Full Name:

interleukin 25

Calculated MW:

177 aa, 20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IP

Species Specificity:

human

Notable Publications

Author	Pubmed ID	Journal	Application
Shu-Yi Yin	27089063	Nat Commun	IP, WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

