

For Research Use Only

RPRD1B Polyclonal antibody, PBS Only

Catalog Number: 24542-1-PBS

Featured Product



Basic Information

Catalog Number:

24542-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21578

GenBank Accession Number:

BC033629

GeneID (NCBI):

58490

UNIPROT ID:

Q9NQG5

Full Name:

regulation of nuclear pre-mRNA domain containing 1B

Calculated MW:

326 aa, 37 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24542-1-AP (RPRD1B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24542-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using 24542-1-AP (RPRD1B antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 24542-1-PBS in a different storage buffer formulation.