

For Research Use Only

TBCEL Polyclonal antibody

Catalog Number: 24548-1-AP



Basic Information

Catalog Number:

24548-1-AP

Size:

150ul, Concentration: 500 ug/ml by Nanodrop and 373 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21725

GenBank Accession Number:

BC120988

GeneID (NCBI):

219899

UNIPROT ID:

Q5QJ74

Full Name:

tubulin folding cofactor E-like

Calculated MW:

424 aa, 48 kDa

Observed MW:

40-48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse testis tissue, PC-3 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 24548-1-AP (TBCEL Antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.