

For Research Use Only

# GALNT13 Polyclonal antibody

Catalog Number: 24604-1-AP



## Basic Information

<b>Catalog Number:</b> 24604-1-AP	<b>GenBank Accession Number:</b> BC101032	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 293 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 114805	
<b>Source:</b> Rabbit	<b>Full Name:</b> UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylglucosaminyltransferase 13 (GalNAc-T13)	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 561 aa, 65 kDa	
<b>Immunogen Catalog Number:</b> AG16526		

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Aliquoting is unnecessary for -20°C storage**

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

