

For Research Use Only

TMPRSS5 Polyclonal antibody

Catalog Number: 24616-1-AP



Basic Information

Catalog Number: 24616-1-AP	GenBank Accession Number: BC121803	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 400 µg/ml by Nanodrop and 267 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 80975	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: transmembrane protease, serine 5	
Isotype: IgG	Calculated MW: 457 aa, 50 kDa	
Immunogen Catalog Number: AG18493	Observed MW: 53 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : rat brain tissue,
Species Specificity: human, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat brain tissue were subjected to SDS PAGE followed by western blot with 24616-1-AP (TMPRSS5 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.