

For Research Use Only

B3GAT2 Polyclonal antibody

Catalog Number: 24631-1-AP



Basic Information

Catalog Number: 24631-1-AP	GenBank Accession Number: BC113995	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 400 µg/ml by Nanodrop;	GeneID (NCBI): 135152	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: beta-1,3-glucuronyltransferase 2 (glucuronosyltransferase S)	
Isotype: IgG	Calculated MW: 323 aa, 37 kDa	
Immunogen Catalog Number: AG18402	Observed MW: 37 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, mouse cerebellum tissue, mouse kidney tissue, mouse retina tissue
Species Specificity: mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24631-1-AP (B3GAT2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.