

For Research Use Only

# GKN2 Polyclonal antibody

Catalog Number: 24651-1-AP



## Basic Information

### Catalog Number:

24651-1-AP

### Size:

150ul , Concentration: 1200 ug/ml by Nanodrop and 600 ug/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG19153

### GenBank Accession Number:

BC146484

### GeneID (NCBI):

200504

### UNIPROT ID:

Q86XP6

### Full Name:

gastrokeine 2

### Calculated MW:

184 aa, 20 kDa

### Observed MW:

20 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

WB 1:500-1:2000

## Applications

### Tested Applications:

WB, ELISA

### Species Specificity:

human, mouse

### Positive Controls:

WB : mouse stomach tissue,

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse stomach tissue were subjected to SDS PAGE followed by western blot with 24651-1-AP (GKN2 Antibody) at dilution of 1:1000 incubated at 4 degree celsius over night.