

For Research Use Only

# TPOR/MPL Polyclonal antibody, PBS Only

Catalog Number: 24760-1-PBS



## Basic Information

Catalog Number:

24760-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19035

GenBank Accession Number:

BC153092

GeneID (NCBI):

4352

UNIPROT ID:

P40238

Full Name:

myeloproliferative leukemia virus oncogene

Calculated MW:

635 aa, 71 kDa

Observed MW:

84 kDa, 74 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Human peripheral blood platelets were subjected to SDS PAGE followed by western blot with 24760-1-AP (TPOR/MPL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24760-1-PBS in a different storage buffer formulation.