

For Research Use Only

SCC4 Polyclonal antibody

Catalog Number: 24770-1-AP



Basic Information

Catalog Number:

24770-1-AP

Size:

150ul , Concentration: 600 ug/ml by Nanodrop and 347 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19943

GenBank Accession Number:

BC063863

GeneID (NCBI):

23383

UNIPROT ID:

Q9Y6X3

Full Name:

KIAA0892

Calculated MW:

613 aa, 69 kDa

Observed MW:

66-69 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PC-3 cells, HeLa cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 24770-1-AP (SCC4 antibody at dilution of 1:600 incubated at room temperature for 1.5 hours.