

For Research Use Only

SLC36A1 Polyclonal antibody, PBS Only

Catalog Number: 24775-1-PBS

Featured Product



Basic Information

Catalog Number:

24775-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20594

GenBank Accession Number:

BC136437

GeneID (NCBI):

206358

UNIPROT ID:

Q7Z2H8

Full Name:

solute carrier family 36 (proton/amino acid symporter), member 1

Calculated MW:

476 aa, 53 kDa

Observed MW:

43-53 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human small intestine tissue slide using 24775-1-AP (SLC36A1 antibody) at dilution of 1:200 (under 0x lens). This data was developed using the same antibody clone with 24775-1-PBS in a different storage buffer formulation.



HeLa cells were subjected to SDS PAGE followed by western blot with 24775-1-AP (SLC36A1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24775-1-PBS in a different storage buffer formulation.