

For Research Use Only

LPPR4 Polyclonal antibody

Catalog Number: 24825-1-AP



Basic Information

Catalog Number:

24825-1-AP

Size:

150ul , Concentration: 700 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20298

GenBank Accession Number:

BC042963

GeneID (NCBI):

9890

UNIPROT ID:

Q7Z2D5

Full Name:

plasticity related gene 1

Calculated MW:

763 aa, 83 kDa

Observed MW:

83-100 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : C6 cells, mouse liver tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24825-1-AP (LPPR4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.