

For Research Use Only

LPPR4 Polyclonal antibody, PBS Only

Catalog Number: 24825-1-PBS



Basic Information

Catalog Number:

24825-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20298

GenBank Accession Number:

BC042963

GeneID (NCBI):

9890

UNIPROT ID:

Q7Z2D5

Full Name:

plasticity related gene 1

Calculated MW:

763 aa, 83 kDa

Observed MW:

83-100 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24825-1-AP (LPPR4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24825-1-PBS in a different storage buffer formulation.