

For Research Use Only

IDNK Polyclonal antibody, PBS Only

Catalog Number: 24831-1-PBS



Basic Information

Catalog Number:

24831-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20088

GenBank Accession Number:

BC142991

GeneID (NCBI):

414328

UNIPROT ID:

Q5T6J7

Full Name:

chromosome 9 open reading frame 103

Calculated MW:

230 aa, 25 kDa

Observed MW:

21 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

The IDNK gene encodes human gluconokinase, which is involved in the phosphorylation of gluconate.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24831-1-AP (IDNK antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24831-1-PBS in a different storage buffer formulation.