

For Research Use Only

LRRC4 Polyclonal antibody, PBS Only

Catalog Number: 24834-1-PBS



Basic Information

Catalog Number:

24834-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20376

GenBank Accession Number:

BC111561

GeneID (NCBI):

64101

UNIPROT ID:

Q9HBW1

Full Name:

leucine rich repeat containing 4

Calculated MW:

653 aa, 73 kDa

Observed MW:

73 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24834-1-AP (LRRC4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24834-1-PBS in a different storage buffer formulation.