

For Research Use Only

SZRD1 Polyclonal antibody

Catalog Number: 24844-1-AP



Basic Information

Catalog Number: 24844-1-AP	GenBank Accession Number: BC023988	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 600 µg/ml by Nanodrop and 287 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 26099	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: chromosome 1 open reading frame 144	
Isotype: IgG	Calculated MW: 152 aa, 17 kDa	
Immunogen Catalog Number: AG20507	Observed MW: 18-22 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2 cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24844-1-AP (SZRD1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.