

For Research Use Only

DTWD2 Polyclonal antibody

Catalog Number: 24873-1-AP



Basic Information

Catalog Number:
24873-1-AP

Size:
150ul , Concentration: 400 µg/ml by Nanodrop and 320 µg/ml by Bradford method using BSA as the standard;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG21665

GenBank Accession Number:
BC114955

GeneID (NCBI):
285605

Full Name:
DTW domain containing 2

Calculated MW:
298 aa, 33 kDa

Observed MW:
33 kDa

Purification Method:
Antigen affinity purification

Recommended Dilutions:
IP 0.5-4.0 ug for IP and 1:200-1:1000 for WB

Applications

Tested Applications:
IP, ELISA

Species Specificity:
human

Positive Controls:

IP : COLO 320 cells,

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP Result of anti-DTWD2 (IP:24873-1-AP, 3ug;
Detection:24873-1-AP 1:300) with COLO 320 cells
lysate 2000ug.