

For Research Use Only

DTWD2 Polyclonal antibody

Catalog Number: 24873-1-AP

1 Publications



Basic Information

Catalog Number:

24873-1-AP

Size:

150ul, Concentration: 400 ug/ml by Nanodrop and 320 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21665

GenBank Accession Number:

BC114955

GeneID (NCBI):

285605

UNIPROT ID:

Q8NBA8

Full Name:

DTW domain containing 2

Calculated MW:

298 aa, 33 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

IP: COLO 320 cells,

Notable Publications

Author	Pubmed ID	Journal	Application
Yixuan Xie	39173631	Cell	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-DTWD2 (IP:24873-1-AP, 3ug;
Detection:24873-1-AP 1:300) with COLO 320 cells
lysate 2000ug.