

For Research Use Only

C1orf161 Polyclonal antibody

Catalog Number: 24886-1-AP



Basic Information

Catalog Number: 24886-1-AP	GenBank Accession Number: BC128148	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 700 µg/ml by Nanodrop;	GeneID (NCBI): 126868	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: chromosome 1 open reading frame 161	
Isotype: IgG	Calculated MW: 362 aa, 42 kDa	
Immunogen Catalog Number: AG21677	Observed MW: 42-45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse skin tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24886-1-AP (C1orf161 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.