

For Research Use Only

TMEM190 Polyclonal antibody

Catalog Number: 24904-1-AP



Basic Information

Catalog Number: 24904-1-AP	GenBank Accession Number: BC128189	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 µg/ml by Nanodrop and 220 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 147744	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	Full Name: transmembrane protein 190	
Isotype: IgG	Calculated MW: 177 aa, 19 kDa	
Immunogen Catalog Number: AG19131		

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human colon cancer tissue, human liver cancer tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 24904-1-AP (TMEM190 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 24904-1-AP (TMEM190 Antibody) at dilution of 1:50 (under 40x lens).