

For Research Use Only

FAM82B Polyclonal antibody, PBS Only

Catalog Number: 24919-1-PBS



Basic Information

Catalog Number:

24919-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20679

GenBank Accession Number:

BC009671

GeneID (NCBI):

51115

UNIPROT ID:

Q96DB5

Full Name:

family with sequence similarity 82, member B

Calculated MW:

314 aa, 36 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

FAM82B, also known as RMDN1 or CGI90, belongs to the FAM82/RMD family. FAM82B is a microtubule-associated protein with potential roles in regulating microtubule dynamics and other cellular functions.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24919-1-AP (FAM82B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24919-1-PBS in a different storage buffer formulation.

IP result of anti-FAM82B (IP:24919-1-AP, 4ug; Detection:24919-1-AP 1:600) with HeLa cells lysate 1760 ug. This data was developed using the same antibody clone with 24919-1-PBS in a different storage buffer formulation.