

For Research Use Only

MCUR1 Polyclonal antibody, PBS Only

Catalog Number: 24948-1-PBS



Basic Information

Catalog Number:

24948-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20040

GenBank Accession Number:

BC016850

GeneID (NCBI):

63933

UNIPROT ID:

Q96AQ8

Full Name:

coiled-coil domain containing 90A

Calculated MW:

359 aa, 40 kDa

Observed MW:

40-50 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24948-1-AP (MCUR1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24948-1-PBS in a different storage buffer formulation.



IP result of anti-MCUR1 (IP:24948-1-AP, 3ug; Detection:24948-1-AP 1:300) with HEK-293 cells lysate 920ug. This data was developed using the same antibody clone with 24948-1-PBS in a different storage buffer formulation.