

For Research Use Only

SPATA1 Polyclonal antibody

Catalog Number: 24980-1-AP



Basic Information

Catalog Number: 24980-1-AP	GenBank Accession Number: BC127634	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 600 µg/ml by Nanodrop and 280 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 64173	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: spermatogenesis associated 1	
Isotype: IgG	Calculated MW: 437 aa, 50 kDa	
Immunogen Catalog Number: AG21784	Observed MW: 40-50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human testis tissue, mouse testis tissue
Species Specificity: human, mouse	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 24980-1-AP (SPATA1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.