

For Research Use Only

AKAP4 Polyclonal antibody, PBS Only

Catalog Number:24986-1-PBS



Basic Information

Catalog Number: 24986-1-PBS	GenBank Accession Number: BC126252	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 8852	
Source: Rabbit	UNIPROT ID: Q5JQC9	
Isotype: IgG	Full Name: A kinase (PRKA) anchor protein 4	
Immunogen Catalog Number: AG21762	Calculated MW: 854 aa, 94 kDa	
	Observed MW: 82 kDa	

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human

Background Information

AKAP4, also named as A-kinase anchor protein 82 kDa, is a 854 amino acid protein, which belongs to the AKAP110 family. AKAP4 localizes to the principle piece of the sperm flagellum and is only expressed in round spermatids. AKAP4 is a major structural component of sperm fibrous sheath and plays a role in sperm motility.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 24986-1-AP (AKAP4 Antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24986-1-PBS in a different storage buffer formulation.