

For Research Use Only

Dexas1 Polyclonal antibody, PBS Only

Catalog Number:25010-1-PBS



Basic Information

Catalog Number:

25010-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18818

GenBank Accession Number:

BC018041

GeneID (NCBI):

51655

UNIPROT ID:

Q9Y272

Full Name:

RAS, dexamethasone-induced 1

Calculated MW:

281 aa, 32 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data