

For Research Use Only

# FAM46C Polyclonal antibody

Catalog Number: 25038-1-AP

Featured Product

2 Publications



## Basic Information

### Catalog Number:

25038-1-AP

### Size:

150ul, Concentration: 400 ug/ml by Nanodrop and 300 ug/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG19147

### GenBank Accession Number:

BC131726

### GeneID (NCBI):

54855

### UNIPROT ID:

Q5VWP2

### Full Name:

family with sequence similarity 46, member C

### Calculated MW:

391 aa, 45 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Cited Applications:

WB

### Species Specificity:

human

### Cited Species:

human

## Notable Publications

| Author        | Pubmed ID | Journal    | Application |
|---------------|-----------|------------|-------------|
| Yuan Xiao Zhu | 28619709  | Cancer Res | WB          |
| Jo Kanasugi   | 32176823  | Cancer Sci | WB          |

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

