

For Research Use Only

ARRDC3 Polyclonal antibody, PBS Only

Catalog Number: 25060-1-PBS

Featured Product



Basic Information

Catalog Number:

25060-1-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18717

GenBank Accession Number:

BC053619

GeneID (NCBI):

57561

UNIPROT ID:

Q96B67

Full Name:

arrestin domain containing 3

Calculated MW:

414 aa, 46 kDa

Observed MW:

46 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data