

For Research Use Only

ARSI Polyclonal antibody, PBS Only

Catalog Number: 25112-1-PBS



Basic Information

Catalog Number:

25112-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18884

GenBank Accession Number:

BC129995

GeneID (NCBI):

340075

UNIPROT ID:

Q5FYB1

Full Name:

arylsulfatase family, member I

Calculated MW:

569 aa, 64 kDa

Observed MW:

64 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 25112-1-AP (ARSI antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25112-1-PBS in a different storage buffer formulation.