

For Research Use Only

AWAT2 Polyclonal antibody

Catalog Number: 25119-1-AP



Basic Information

Catalog Number: 25119-1-AP	GenBank Accession Number: BC153070	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 240 µg/ml by Nanodrop and 207 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 158835	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: acyl-CoA wax alcohol acyltransferase 2	
Isotype: IgG	Calculated MW: 333 aa, 38 kDa	
Immunogen Catalog Number: AG18936	Observed MW: 38-42 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 25119-1-AP (AWAT2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.