

For Research Use Only

acetylated Tubulin (Lys40) Polyclonal antibody



Catalog Number: 25126-1-AP

2 Publications

Basic Information

Catalog Number: 25126-1-AP	GenBank Accession Number: NM_006009	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 193 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 7846	
Source: Rabbit	Full Name: tubulin, alpha 1a	
Isotype: IgG	Calculated MW: 52 kDa	
	Observed MW: 50-55 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IF, WB

Species Specificity:
human, mouse, rat

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Hidetoshi Tanioka	35861932	Jpn J Ophthalmol	IF
Luojia Liu	37949303	Biochim Biophys Acta Mol Cell Res	WB,IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data