

For Research Use Only

# LAMB4 Polyclonal antibody

Catalog Number: 25217-1-AP



## Basic Information

### Catalog Number:

25217-1-AP

### Size:

150ul , Concentration: 900 ug/ml by Nanodrop and 380 ug/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG19401

### GenBank Accession Number:

BC140804

### GeneID (NCBI):

22798

### UNIPROT ID:

A4D0S4

### Full Name:

laminin, beta 4

### Calculated MW:

1761 aa, 194 kDa

### Observed MW:

193 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

WB 1:500-1:1000

## Applications

### Tested Applications:

WB, ELISA

### Species Specificity:

human, rat, mouse

### Positive Controls:

WB : mouse skin tissue,

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 25217-1-AP (LAMB4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.