

For Research Use Only

DACH2 Polyclonal antibody, PBS Only

Catalog Number: 25219-1-PBS



Basic Information

Catalog Number:

25219-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18403

GenBank Accession Number:

BC114950

GeneID (NCBI):

117154

UNIPROT ID:

Q96NX9

Full Name:

dachshund homolog 2 (Drosophila)

Calculated MW:

599 aa, 65 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Dachshund homolog 2 (DACH2), is a 599 amino acid protein, which belongs to the DACH family. DACH2 localizes in the nucleus and as a transcription factor is involved in regulation of organogenesis.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25219-1-AP (DACH2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25219-1-PBS in a different storage buffer formulation.