

For Research Use Only

IRGM Polyclonal antibody, PBS Only

Catalog Number: 25231-1-PBS



Basic Information

Catalog Number:

25231-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18631

GenBank Accession Number:

BC128168

GeneID (NCBI):

345611

UNIPROT ID:

A1A4Y4

Full Name:

immunity-related GTPase family, M

Calculated MW:

181 aa, 20 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25231-1-AP (IRGM antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25231-1-PBS in a different storage buffer formulation.