

For Research Use Only

NTAQ1/WDYHV1 Polyclonal antibody

Catalog Number: 25308-1-AP



Basic Information

Catalog Number:

25308-1-AP

Size:

150ul, Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14713

GenBank Accession Number:

BC008781

GeneID (NCBI):

55093

UNIPROT ID:

Q96HA8

Full Name:

WDYHV motif containing 1

Calculated MW:

205 aa, 24 kDa

Observed MW:

24 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : MCF-7 cells, U-251 cells, U2OS cells, mouse pancreas tissue, rat pancreas tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25308-1-AP (WDYHV1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.